



HD 3406.2 BENCH-TOP CONDUCTIVITY METER

The HD3406.2 is a bench top instrument for electrochemical measures: **conductivity and temperature**.

The displayed data can be stored (**datalogger**) and can be transferred to PC or serial printer thanks to the multi-standard serial port RS232C and USB2.0 and software DeltaLog9 (Vers.2.0 and subsequent ones). The storing and printing parameters can be set from menu.

The HD3406.2 measures **conductivity, liquid resistivity in liquids, total dissolved solids (TDS) and salinity** using combined 4-ring and 2-ring conductivity/temperature probes. Temperature is measured by Pt100 or Pt1000 immersion, penetration or contact probes.

The probe calibration can be performed automatically in one or more of the 147 μ S, 1413 μ S, 12880 μ S or 111800 μ S/cm conductivity calibration solutions.

The display shows continually the temperature in °C or °F and one selectable parameter according to the connected probe type, i.e. in case of conductivity probe it is possible to select between χ or Ω or TDS or NaCl.

Other functions of this instrument include: Max, Min and Avg function, the Auto-HOLD function, the automatic turning off which can also be excluded.

The instruments have IP66 protection degree.



Technical characteristics HD3406.2 χ , Ω , TDS, NaCl, °C/°F measurement

Instrument

Dimensions (Length x Width x Height)

Weight

Materials

Display

220x120x55mm

460g (complete with batteries)

ABS, rubber

2x4½ characters plus symbols

visible area: 52x42mm

Operating conditions

Working temperature

Storage temperature

Working relative humidity

Protection degree

-5 ... 50°C

-25 ... 65°C

0 ... 90% RH without condensation

IP66

Power

Batteries

Autonomy (only batteries)

Mains (cod. SWD10)

3 batteries 1.5V type AA

100 hours with 1800mAh alkaline batteries

Output mains adapter 100-240Vac/ 12Vdc-1A

Security of memorized data

Unlimited

Selectable storage interval

1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

Time

Date and hour

Accuracy

Schedule in real time

1min/month max drift

Serial interface RS232C

Type

Baud rate

Data bit

Parity

Stop bit

Flow Control

Serial cable length

Selectable print interval

RS232C electrically isolated

Can be set from 1200 to 38400 baud

8

None

1

Xon/Xoff

Max 15m

immediate or 1s, 5s, 10s, 15s, 30s, 1min, 2min, 5min, 10min, 15min, 20min, 30min and 1hour

USB Interface

Type

1.1 - 2.0 electrically isolated

Common connections to all models

Serial interface and USB

Mains adapter (cod. SWD10)

8-pole MiniDin connector

2-pole connector (positive at centre) 12Vdc/1A

Storage of measured values

Type

Quantity

2000 pages of 18 samples each

36,000 sets of measures made up of [χ - Ω or TDS or NaCl] and [°C - °F]

Measurement connections

Input conductivity

Input for temperature probes

complete with TP47 modules

8-pole male DIN45326 connector

8-pole male DIN45326 connector

Measurement of conductivity by instrument

Measurement range (Kcell=0.01)

Measurement range (Kcell=0.1)

Measurement range (Kcell=1)

Resolution

0.000...1.999 μ S/cm 0.001 μ S/cm

0.00...19.99 μ S/cm 0.01 μ S/cm

0.0...199.9 μ S/cm 0.1 μ S/cm

200...1999 μ S/cm 1 μ S/cm

2.00...19.99mS/cm 0.01mS/cm


20.00...199.9mS/cm 0.1mS/cm

Vdc in

N.C.

8-pole Mini DIN

8-pole DIN 45326



Range di misura (Kcell=10).	200...1999mS/cm	1mS/cm
Accuracy (conductivity)	$\pm 0.5\% \pm 1\text{digit}$	
Measurement of resistivity by instrument		
Measurement range (Kcell=0.01)	Up to 1G Ω ·cm	(*)
Measurement range (Kcell=0.1)	Up to 100M Ω ·cm	(*)
Measurement range (Kcell=1)	5.0...199.9 Ω ·cm	0.1 Ω ·cm
	200...999 Ω ·cm	1 Ω ·cm
	1.00k...19.99k Ω ·cm	0.01k Ω ·cm
	20.0k...99.9k Ω ·cm	0.1k Ω ·cm
	100k...999k Ω ·cm	1k Ω ·cm
	1...10M Ω ·cm	1M Ω ·cm
Measurement range (Kcell=10)	0.5...5.0 Ω ·cm	0.1 Ω ·cm
Accuracy (resistivity)	$\pm 0.5\% \pm 1\text{digit}$	
Measurement of total dissolved solids (with coefficient $\chi/TDS=0.5$)		
Measurement range (Kcell=0.01)	0.00...1.999mg/l	0.005mg/l
Measurement range (Kcell=0.1)	0.00...19.99mg/l	0.05mg/l
Measurement range (Kcell=1)	0.0...199.9 mg/l	0.5 mg/l
	200...1999 mg/l	1 mg/l
	2.00...19.99 g/l	0.01 g/l
	20.0...99.9 g/l	0.1 g/l
Measurement range (Kcell=10)	100...999 g/l	1 g/l
Accuracy (total dissolved solids)	$\pm 0.5\% \pm 1\text{digit}$	
Measurement of salinity		
Measurement range / Resolution	0.000...1.999g/l	1mg/l
	2.00...19.99g/l	10mg/l
	20.0...199.9g/l	0.1g/l
Accuracy (salinity)	$\pm 0.5\% \pm 1\text{digit}$	
Temperature measurement by instrument		
Measurement range Pt100	-50...+200°C	
Measurement range Pt1000	-50...+200°C	
Resolution	0.1°C	
Accuracy	$\pm 0.25^\circ\text{C}$	
Drift after 1 year	0.1°C/year	
Automatic/manual temperature compensation		
Reference temperature	0...100°C with $\alpha_r = 0.00...4.00\%/\text{°C}$	
Conversion factor χ/TDS	20°C or 25°C selectable from menu	
Cell constant K (cm^{-1})	0.4...0.8	
	0.01 - 0.1 - 0.7 - 1.0 - 10.0	
Standard solutions automatically detected (@25°C)		
	147 $\mu\text{S}/\text{cm}$	
	1413 $\mu\text{S}/\text{cm}$	
	12880 $\mu\text{S}/\text{cm}$	
	111800 $\mu\text{S}/\text{cm}$	

(*) The resistivity measurement is obtained from the reciprocal of conductivity measurement. Close to the bottom of the scale, the indication of resistivity appears like reported in the table below:

K cell = 0.01 cm^{-1}		K cell = 0.1 cm^{-1}	
Conductivity ($\mu\text{S}/\text{cm}$)	Resistivity (M Ω ·cm)	Conductivity ($\mu\text{S}/\text{cm}$)	Resistivity (M Ω ·cm)
0.001 $\mu\text{S}/\text{cm}$	1000 M Ω ·cm	0.01 $\mu\text{S}/\text{cm}$	100 M Ω ·cm
0.002 $\mu\text{S}/\text{cm}$	500 M Ω ·cm	0.02 $\mu\text{S}/\text{cm}$	50 M Ω ·cm
0.003 $\mu\text{S}/\text{cm}$	333 M Ω ·cm	0.03 $\mu\text{S}/\text{cm}$	33 M Ω ·cm
0.004 $\mu\text{S}/\text{cm}$	250 M Ω ·cm	0.04 $\mu\text{S}/\text{cm}$	25 M Ω ·cm
...

ORDERING CODES

HD3406.2: The kit is composed of: instrument HD3406.2 **datalogger**, for measurement of conductivity - resistivity - TDS - salinity - temperature, 3 1.5V alkaline batteries, operating manual and **DeltaLog9 version 2.0**.
pH/mV electrodes, conductivity probes, dissolved oxygen probes, temperature probes, standard reference solutions for different measurement types, connection cables for pH electrodes with S7 connector, cables for data download to PC or printer have to be ordered separately.



χ



Ω



mg/l



NaCl

ACCESSORIES

HD2110CSNM: 8-pole connection cable Mini Din - Sub D 9-pole female for RS232C, for connection to PC without USB input.
HD2101/USB: Connection cable USB 2.0 connector type A - 8-pole Mini Din for connection to PC with USB input.
SWD10: Stabilized power supply at 100-240Vac/12Vdc-1A mains voltage.
HD40.1: Portable, serial input, 24 column thermal printer, 57mm paper width.
HD22.2: Laboratory electrode holder composed of basis plate with incorporated magnetic stirrer, staff and replaceable electrode holder. Height max. 380mm.
HD22.3: Laboratory electrode holder with metal basis plate. Flexible electrode holder for free positioning. For Ø 12mm probes.
TP47: Module for the connection of Pt100 4-wire and Pt1000 2-wire probes.

Combined conductivity and temperature probes

SP06T: Combined conductivity and temperature 4-electrode cell in Platinum, body in Pocan. Cell constant K = 0.7. Measurement range 5 $\mu\text{S}/\text{cm}$... 200mS/cm, 0...90°C.

SPT401.001: Combined conductivity and temperature 2-electrode cell in stainless steel AISI 316. Cell constant K = 0.01. Measurement range 0.04 $\mu\text{S}/\text{cm}$... 20 $\mu\text{S}/\text{cm}$, 0...120°C. Measurement in closed-cell..

SPT01G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 0.1. Measurement range 0.1 $\mu\text{S}/\text{cm}$... 500 $\mu\text{S}/\text{cm}$, 0...80°C.

SPT1G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 1. Measurement range 10 $\mu\text{S}/\text{cm}$... 10mS/cm, 0...80°C.

SPT10G: Combined conductivity and temperature 2-electrode Platinum-wire cell, body in glass. Cell constant K = 10. Measurement range 500 $\mu\text{S}/\text{cm}$... 200mS/cm, 0...80°C.

Electrode characteristics at page 402

Standard conductivity calibration solutions

HD8747: Standard calibration solution 0.001mol/l equal to 147 $\mu\text{S}/\text{cm}$ @25°C - 200cc.

HD8714: Standard calibration solution 0.01mol/l equal to 1413 $\mu\text{S}/\text{cm}$ @25°C - 200cc.

HD8712: Standard calibration solution 0.1mol/l equal to 12880 $\mu\text{S}/\text{cm}$ @25°C - 200cc.

HD87111: Standard calibration solution 1mol/l equal to 111800 $\mu\text{S}/\text{cm}$ @25°C - 200cc.

Temperature probes complete with TP47 module

TP47.100: Direct 4 wires Pt100 sensor immersion probe. Probe's stem Ø 3mm, length 230mm. Connection cable 4 wires with connector, length 2 m.

TP47.1000: Pt1000 sensor immersion probe. Probe's stem Ø 3mm, length 230mm. Connection cable 2 wires with connector, length 2 m.

TP87.100: Pt100 sensor immersion probe. Probe's stem Ø 3mm, length 70mm. 4 wire connection cable with connector, length 1 m.

TP87.1000: Pt1000 sensor immersion probe. Probe's stem Ø 3mm, length 70mm. 2 wire connection cable with connector, length 1 m.

Accessories

TP47: Module for the connection of Pt100 4-wire and Pt1000 2-wire probes.